



MINOS Status Report



- Detector and Magnet ON
- Electronics Calibrations issues, lately seen with the new DAQ, have been fixed (Bill Badgett).



MINOS Status Report



- Still investigating some readout issues
 - anomalous erratic Pedestal patterns
 - similar effects in Singles rates plots
- General work plan is to address and fix hardware issues to keep detector ready for beam
 - 3 problems channels (hot/warm/dead)
 - need to replace Minder boards
- Spares situation
 - 6 Minders repaired at Argonne
 - 12 new being assembled at Argonne
 - 1 working VTM : 2 to be repaired (rarely break)



MINOS Status Report



- Near Detector New Networking Rack Update
 - passed ORC review and is powered ON
 - switch over from old rack planned for this Thursday, Jun 27



MINOS Status Report



- Developed readout and DAQ issue on Branch 2 this past Saturday
 - working on resolving it (swap in spare boards etc.)
- RPS alarms on Crate 15
 - working on it



MINOS Status Report



- Regular Shifts are underway